

## Organism Directory

# *Clostridium perfringens*

### Introduction

*Clostridium perfringens* is the most common *Clostridium* found in cases of gas gangrene (clostridial myonecrosis) with a high mortality rate. However, the recovery of this organism from wounds is not necessarily an indication of this condition.

Gas gangrene is a clinical diagnosis and the recovery of this organism merely supports the diagnosis. It causes a wide spectrum of other foodborne illnesses and acute gastroenteritis.

### Appearance

Large, translucent grey colonies surrounded by a double zone of haemolysis on anaerobic blood agar. Colonies may be round or take on a rhizoid appearance. The inner zone of haemolysis is clearer than the outer zone giving a target-like appearance. This double zone is much clearer on sheep blood agar, although this is rarely used in the UK



*C. perfringens* at 48 hours

### Features

- Gram-positive rod bacterium
- Obligate anaerobe
- Spore-former (under specialised conditions)
- Toxin-former
- 35-37°C for 48-72 hours
- Causes gas gangrene
- Catalase negative
- Oxidase negative
- Grows on fastidious anaerobe agar and Bruella blood agar



## Practical & Clinical Microbiology of Anaerobes

Held in Cardiff and delivered by the UK Anaerobe Reference Unit, this 2-day residential course offers the chance to learn from the experts about the latest technologies and techniques. CPD accredited for 11 credits. Invaluable preparation for FRCPath.

For a more details please email:  
Bethany Rooke – [sales@dwscientific.co.uk](mailto:sales@dwscientific.co.uk)



## UKAS Calibration/Validation

DWS is accredited by UKAS to BS EN ISO/IEC 17025:2017 for:

- Temperature calibration within the range -80°C to 140°C (Autoclaves, media preparators, and incubators)
- Non-automatic weighing instrument calibration within the range 0-9kg (Gravimetric diluters and balances)

Email Sally Baines for more information  
– [service@dwscientific.co.uk](mailto:service@dwscientific.co.uk)



## Contract Laboratory

The work of our microbiology laboratories extends from routine analysis/quality control to substantial research and development programmes for major multinational pharmaceutical companies. In addition to our membership of the UK GLP compliance program, we comply with national and international standards as required, including BS / EN / ISO / ASTM test methods, BSAC, EUCAST and CLSI guidelines.

Email Dr Andrew Pridmore to discuss your requirements –  
[laboratory@dwscientific.co.uk](mailto:laboratory@dwscientific.co.uk)



## Exploring the Microbiome

Each installation is only 30 minutes long and the recordings are uploaded to our website!

Visit [www.dwscientific.com/dwswebinars](http://www.dwscientific.com/dwswebinars)  
to begin the exploration!

