

Canopy Technical Specifications

The canopy is RoHS, REACH and FCC compliant, CE certified, made in the USA, and has a two-year warranty.



General

Compatible Products	▶ Alto Microplate Reader
Dimensions	▶ 11.32 X 6.61 X 4.78 cm (LWH)
Weight	▶ 170 g
Power Input	▶ 5.1 V, 1 A via USB-C receptacle

Communication

Wireless Protocol	▶ IEEE 802.15.4
Wireless Range	▶ 15 m

Cerillo Labrador Premium Software

Operating System	▶ macOS: 12 (Monterey); Windows: 10, 11
------------------	---