

IMPORTANT INFORMATION TO READ and RETURN

Installation Requirements for DILUWEL

Thank you for choosing one of our products for your laboratory. To enable our engineers to perform an efficient, trouble-free installation please study, complete and email this form to us on **service@dwscientific.co.uk**. Should you have any questions, please do not hesitate to contact us, as we are here to help. When we have received the completed form, our Service Department will contact you to arrange a mutually convenient installation date.

The following information represents the ideal requirement.

Please contact us IMMEDIATELY if your intended location does not match this specification.

Tica	se comact us	SIMINIEDIATEE	in your interded location does not match this specification	✓
Space Requir	ements			
The weight of th	e equipment is 8	kg. The bench allo	cated must be flat, level and of sufficient size to support the base fully.	
External Dim	ensions			
Width	Depth	Height		
mm	mm	mm		
250	275	400		
Environments Altitude	al Conditions	up to 2,000m from 0 to 45°		
Temperature Relative humidity Machine supply voltage Power rating		from 10 to 45 from 10 to 80 24 V DC, 3.7 120W	%	
Mains Requir	ements			
Electricity Supply		Socket		
100 – 230 CAV, 50-60Hz		1 x 3 pin, 13 ar	np	

Delivery and installation are free of charge (unless otherwise agreed). If our engineers are unable to install the unit and a return journey is necessary, **a charge may be made**.

It is essential that this form is completed and returned to avoid delay to your installation.

THANK YOU FOR THINKING WHITLEY

Signature	Title	
Print Name	Establishment	

DWS/SSP0005 Form Revision Level: 0