

IMPORTANT INFORMATION TO READ and RETURN

Installation Requirements for the DISTRIWEL 440

Thank you for choosing one of our products for your laboratory. To enable our engineers to perform an efficient, trouble-free installation please study, complete and email this form to us on **service@dwscientific.co.uk**. Should you have any questions, please do not hesitate to contact us, as we are here to help. When we have received the completed form, our Service Department will contact you to arrange a mutually convenient installation date.

**The following information represents the ideal requirement.
Please contact us IMMEDIATELY if your intended location does not match this specification.**



Access Requirements	For access, the dimensions stated below should be taken into account when checking the size of lifts, stairwells, doorways, etc.	<input type="checkbox"/>									
Space Requirements	The weight of the equipment is 30kg. The bench allocated must be flat, level and of sufficient size to support the base fully. The unit should be sited as close as possible (preferably to the left) of the Media Preparator.										
Bench Space Required	<table border="1" style="width: 100%; border-collapse: collapse; text-align: center;"> <thead> <tr> <th style="width: 33%;">Width</th> <th style="width: 33%;">Depth</th> <th style="width: 33%;">Height</th> </tr> <tr> <td><i>mm</i></td> <td><i>mm</i></td> <td><i>mm</i></td> </tr> <tr> <td>808</td> <td>689</td> <td>940</td> </tr> </thead> </table>	Width	Depth	Height	<i>mm</i>	<i>mm</i>	<i>mm</i>	808	689	940	<input type="checkbox"/>
Width	Depth	Height									
<i>mm</i>	<i>mm</i>	<i>mm</i>									
808	689	940									
Mains Requirements	<table style="width: 100%;"> <tr> <td style="width: 50%;">Electricity Supply</td> <td style="width: 50%;">Socket</td> </tr> <tr> <td>230 volts +/-10%, 50-60Hz or 110 volts +/-10%, 50-60Hz 150W</td> <td>1 x Three pin, 13 amp</td> </tr> </table>	Electricity Supply	Socket	230 volts +/-10%, 50-60Hz or 110 volts +/-10%, 50-60Hz 150W	1 x Three pin, 13 amp	<input type="checkbox"/>					
Electricity Supply	Socket										
230 volts +/-10%, 50-60Hz or 110 volts +/-10%, 50-60Hz 150W	1 x Three pin, 13 amp										

Delivery and installation are free of charge (unless otherwise agreed). If our engineers are unable to install the unit and a return journey is necessary, **a charge may be made.**

It is essential that this form is completed and returned to avoid delay to your installation.

THANK YOU FOR THINKING WHITLEY

Signature	<input style="background-color: white; color: black;" type="text"/>	Title	<input style="background-color: white; color: black;" type="text"/>
Print Name	<input style="background-color: white; color: black;" type="text"/>	Establishment	<input style="background-color: white; color: black;" type="text"/>