

IMPORTANT INFORMATION TO READ and RETURN

Installation Requirements for MIXWEL

Thank you for choosing one of our products for your laboratory. To enable our engineers to perform an efficient, trouble-free installation please study, complete and email this form to us on **service@dwscientific.co.uk**. Should you have any questions, please do not hesitate to contact us, as we are here to help. When we have received the completed form, our Service Department will contact you to arrange a mutually convenient installation date.

**The following information represents the ideal requirement.
Please contact us IMMEDIATELY if your intended location does not match this specification.**



Space Requirements

The weight of the equipment is 14kg. The bench allocated must be flat, level and of sufficient size to support the base fully.

External Dimensions

Width	Depth	Height
<i>mm</i>	<i>mm</i>	<i>mm</i>
220	410	340

Mains Requirements

Electricity Supply

100 - 230 Volts (+/- 10%) 50-60Hz 0.5A

Power Output: 120W

Socket

1 x Three Pin, 13 amp

Delivery and installation are free of charge (unless otherwise agreed). If our engineers are unable to install the unit and a return journey is necessary, **a charge may be made.**

It is essential that this form is completed and returned to avoid delay to your installation.

THANK YOU FOR THINKING WHITLEY

Signature

Title

Print Name

Establishment