

## FACT SHEET

### MASTERCLAVE 09 MEDIA PREPARATOR

The Masterclave 09 can produce up to 9 litres of media in around 60 minutes. It is easy to use and can be linked to plate pouring devices and gravimetric diluters.



#### Features and Benefits

- At full capacity this media preparator can produce 100 x 90ml dilution broths or 530 x 17ml agar plates: larger batches means less QC.
- Produce high quality media due to microprocessor controlled heating (accurate to within  $\pm 1^{\circ}\text{C}$ ) with no danger of overheating.
- Lightweight inner vessel for easy removal and cleaning.
- Eliminates the risks associated with manual preparation of media in glass vessels.
- Optional chart recorder allows full traceability.
- Optional fully comprehensive maintenance and breakdown contracts are available to prolong the life of your investment. (Note: different arrangements apply in different countries).
- Optional calibration and validation for UKAS.

#### Facts and Figures

**Weight:** 60kg

**Dimensions:** 435 x 500 x 590 mm (D x L x H)

**Electricity Requirements:** 230  $\pm$  10% V AC or 115  $\pm$  10% V AC – single phase