

## FACT SHEET

### SMASHER®

The quietest peristaltic laboratory blender on the market, Smasher® produces perfectly homogenised samples.



### Features and Benefits

- Proven to operate at 10dBA less than other standard blenders – **8 times quieter** – so is suitable for use anywhere in your laboratory.
- No jam, patented smash effect to ensure continuous operation.
- Quick release paddles that can be removed without the use of tools and a totally accessible chamber make cleaning an easy task.
- Easy to use, digital keypad for the programming of different blending cycles (from 10 to 180 seconds or continuously).
- LCD display screen allows you to constantly monitor the blending duration and speed.
- Blending capacity of 80ml to 400ml bags giving you plenty of flexibility.
- ISO7518 and ISO6887 compliant.
- Optional fully comprehensive maintenance and breakdown contracts are available to prolong the life of your investment. (Note: different arrangements apply in different countries).

### Facts and Figures

**Weight:** 21 kg

**Dimensions:** 380 x 400 x 330 mm (D x L x H)

**Electricity Requirements:** 230 ±10% V AC or 115 ± 10% V AC - single phase