

NEW



**don whitley
scientific**
www.dwscientific.com

Whitley H25 Hypoxystation

Don Whitley Scientific is proud to launch their most compact cell culture workstation, the Whitley H25 Hypoxystation. The H25 is equipped with many unique features to ensure easy and efficient use and to guarantee that the very best conditions are maintained. It has a capacity of 145 x 96 well plates while still leaving more than adequate space to process samples in a strictly controlled environment. The H25 allows you to control oxygen in 0.1% increments up to 20%; carbon dioxide in 0.1% increments up to 15%; and relative humidity up to 80%.

This workstation is equipped with sleeved ports and has a built-in, 6 litre rapid airlock that can accommodate up to 15 x 96 well plates or 3 x 500 ml Duran bottles. The airlock allows items to be transferred into the chamber in less than 20 seconds and the doors are interlocked so that inner and outer doors cannot be opened simultaneously. A side loading letterbox entry system is also available, providing a straightforward way to quickly introduce individual samples into the workstation. The H25 is fitted with a touch screen providing the user with an intuitive interface to display the parameters and control the workstation.

With a wide selection of options available, the H25 can fulfil the needs of any research laboratory requiring precise control to create physiologically relevant conditions.

Dimensions (w/d/h) - 790mm/31" x 720mm/28.5" x 710mm/28"

Weight - 74kg/163lb