

IMPORTANT INFORMATION – PLEASE READ

Installation Requirements for WASP Spiral Plater

Thank you for choosing one of our products for your laboratory. To enable our engineers to perform an efficient, trouble-free installation please study, complete and fax this form to us on **01274 531197**.

Should you have any questions, please do not hesitate to contact us, as we are here to help.

When we have received the completed form, our Service Department will contact you to arrange a mutually convenient installation date.

Please contact us if you have any further questions about the installation procedure.



1. Do you have a vacuum source or have you ordered the 602 Vacuum Source?	<input type="checkbox"/>
2. Have you allocated flat, level benchspace that can support 17.5kg and an item of 410 x 530 x 360mm (D x L x H)?	<input type="checkbox"/>
3. Have you allocated two spare mains outlets for the unit (230 +/- 10%V AC)?	<input type="checkbox"/>
4. Do you have suitable disinfectant available to fill the bottles?	<input type="checkbox"/>
5. Have you seen Technical Note WA06 which sets out the best way to prepare agar plates for the WASP?	<input type="checkbox"/>
6. Have you ordered sufficient disposable beakers?	<input type="checkbox"/>
7. Are you fully aware of the Quality Counts Scheme and what actions are required of you (UK and Ireland only)?	<input type="checkbox"/>
8. Are you using Petri dishes from more than one supplier? - Have you been made aware of the adjustable turntable option?	Yes / No Yes / No
9. Will real samples be put through the WASP on the day of installation?	Yes / No
10. If you are manually counting plates, is any additional training required?	Yes / No
11. Are you using pre-poured media? If so, who is the supplier?	Yes / No
<input style="width: 50%; height: 20px;" type="text"/>	

In the UK, delivery and installation are free of charge (unless otherwise agreed). If our engineers are unable to install the unit and a return journey is necessary, **a charge may be made**. Export customers, please refer to your local distributor.

It is essential that this form is completed and returned, to avoid delay to your installation.

THANK YOU FOR THINKING WHITLEY

Signature	<input style="width: 95%; height: 20px;" type="text"/>	Title	<input style="width: 95%; height: 20px;" type="text"/>
Print Name	<input style="width: 95%; height: 20px;" type="text"/>	Establishment	<input style="width: 95%; height: 20px;" type="text"/>